



# LiteShower™ Portable Handicap Shower Stalls

**Setup in minutes, take showers anywhere at home**

## How It Works

1. Hand-held shower head with 7' hose connects to indoor standard faucet (quick-connect adapter included);
2. Shower tray sits on floor. Tray has flexible walls that bend flat under wheelchairs but bounce back afterwards;
3. (2) free-standing shower screens snap together around the tray. Shower curtains open and close on all four sides;
4. Waste water in the tray is pumped to sink. Pump sits outside the tray and has suction hose going into the tray.



shower curtain opens and closes on all four sides



LiteShower Recliner Model



LiteShower™ Shower-in-Bed (patent-pending)

## Models

System	UPC Code	Overall Size	Folded Screen Size	Price
LiteShower Standard	719926241878	3' wide x 4' long x 4' tall	4' long x 4' tall x 2" thick each	\$799
LiteShower Recliner	719926241892	3' wide x 7' long x 3' tall	4' long x 3' tall x 3" thick each	\$899
LiteShower Tall	719926241885	3' wide x 4' long x 6' tall	6' long x 4' tall x 2" thick each	\$899
LiteShower Bariatrics	719926241915	4' wide x 4' long x 4' tall	4' long x 4' tall x 2" thick each	\$899
LiteShower Shower-in-Bed	719926241908	3' wide x 7' long, frameless	No screens needed	\$799

**FREE delivery with Fedex Ground**

Complete set includes (screens not needed for Shower-in-Bed)

Assembly Requirement:

- (1) handheld shower head with 7' hose and on/off control
- (1) adaptor to connect hose to faucet
- (1) shower divert valve to connect hose to existing shower
- (1) self-support, no setup, wheelchair-accessible shower tray
- (1) pair of folding screens with sliding curtains and hooks
- (1) pump with hose, filter to dispose waster water into sink

- 1. Replace threaded aerator with faucet adaptor (one-time)
- 2. Assemble screens with a screwdriver (one-time)
- 3. Hand-tight inline filter, fitting, and hose to pump (one-time)

**1-800-561-3576**

[www.LiteShower.com](http://www.LiteShower.com)

**Made in USA**